



Road Five, Winsford Industrial Estate,
Winsford, Cheshire CW7 3RB

Tel: +0044 1606 593456
Fax: +0044 1606 861354

Web: www.renrayhealthcare.com
Email: sales@renrayhealthcare.com



NHS Contract Reference:
2013/S 022-033974

This brochure is printed on a FSC paper produced at a mill that has been awarded the ISO14001 certificate for environmental management.
Printed by a FSC certified printer using vegetable based inks.





About Renray

At Renray we have been producing high quality furniture for over 45 years and are the UK's largest and leading supplier to the healthcare sector.

Whether you require just a fast efficient delivery of quality furniture or a full room installation and fitting service, we have the experience and resources to handle your contract. Projects undertaken for the NHS and most leading independent healthcare groups range from one off units to a seven hundred plus bedroom accommodation centre, all of which receive the same attention to detail and quality service.

We manufacture and assemble in our own purpose built factory in Winsford, Cheshire, as well as in Europe to British Standards. Hence we are able to ensure that your furniture is produced to the highest standards and we work with you to arrange, plan and meet your time schedule and budget. Our own fleet of vehicles and staff deliver and install all of our furniture.

We understand that you are purchasing furniture for the long term and that you require products that are both fit for purpose, stylish and will continue to perform their function well into the future, which is why we design and build our furniture with you in mind.

Contents

About us

About Renray	2
Recognised for quality craftsmanship	4
Consideration for the global environment	5
Total peace of mind due to our business controls	5
A full turnkey solution	6
And just to make things easier for you	7
Specially designed for use by the NHS	8
High performance fabrics for a ward environment	9
Express Delivery	9

Patient bedside furniture

Enterprise tables	12 - 13
Enterprise lockers	14 - 15
600 range lockers	16
A range lockers	17
100 range combination lockers	18
Elite patient chair range	19 - 21
Baroque chair range	22
Balero chair range	23
Visitor seating	24

Bariatric / recliners & specialist seating

Bariatric chair range	26
Recliners	27
Specialist seating	28

Waiting area

Finland chairs (stackable)	30
Venetian chairs (stackable)	31
Venetian chairs	32
Baroque chairs	33
Balero chairs	34
Coffee tables	35

Dining furniture

Bistro range	38 - 39
Chelford dining chairs	40
Dining tables	41

Sofas, upholstered chairs & tub style seating

Belton range	44
Harvard range	45
Milan range	46
Kingsley range	47
Chelford range	48
Oulton range	49

Your next step

Contact us	51
------------------	----

We also have a full range of bedroom furniture and beds - contact us for details

Recognised for quality craftsmanship

At Renray our qualified craftsmen manufacture furniture to meet the high standards we set and, where appropriate, meet British Standards. These nationally recognised standards include:

- BS 5852 Part 2, Crib 5 Flame Retardancy
- BS EN 15373:2007 Strength & Durability
- BS EN 7176:2007 Resistance to Ignition
- BS 5867-2 type B & C Fabric for curtains, drapes & blinds
- BS 7177:2008 Resistance against ignition of mattresses, mattress pads, divans and bed bases

Consideration for the global environment



Our policy is to provide and maintain working practices that have the least detrimental effect on the global environment and therefore at Renray we are FSC accredited. We operate a lean manufacturing process whereby all our products are made to order by our skilled workforce.



Total peace of mind due to our business controls



When ordering any Renray product you have the peace of mind of knowing that we have been audited and accredited and have the necessary standards and business controls in place. We have achieved **ISO 14001** for our environmental management systems, **OHSAS 18001** for our health and safety management systems, and **ISO 9001** for our quality management systems.



A full turnkey solution

If you require it, we are specialists in providing a complete turnkey package ranging from chairs and tables, a full range of patient bedside furniture plus specialised seating, and beds and mattresses, through to carpets and curtains - in fact everything for the ward environment. The package also includes free room placement, free delivery, full project management, relevant product training and removal of waste materials.

And just to make things easier for you...



Visit our showrooms

To help you make your decision, visit our showrooms where you can view and try out our furniture. The showrooms are located at our office, just off junction 18 of the M6 in Winsford, Cheshire. To arrange a visit please contact our Sales Office on **01606 593456**.

Order delivery

Delivery is generally between four and six weeks after receipt of your order. Please enquire if a quicker delivery is required (see page 9 for Express Delivery details). All major contracts will be placed free of charge in the rooms specified by the customer, providing reasonable access is available. All orders are delivered by our own fleet and our own staff, all of whom are CRB checked.

Quotations and samples

For a competitive quotation please contact our Sales Office on **01606 593456**. One of our nationally based Area Sales Managers will be happy to arrange a meeting at a location of your choice to discuss your needs. Furniture samples can be made available to you and delivered for a two week trial period.

After sales service

If you have any further needs after your order has been delivered, or any comments and suggestions, please contact us. We value the ongoing relationship we have with our customers.

NHS Contract Reference:
2013/S 022-033974

Specially designed for use by the NHS

Developed specifically for the NHS, our range of furniture has been thoughtfully designed, carefully constructed and rigorously tested to meet the needs of today's hospitals.

Our range has been created after research with healthcare providers and with their patients in mind. We have worked to create safe, robust and attractive furniture which can form part of a modern, state of the art hospital environment.

High performance fabrics for a ward environment

All of our chairs featured in this brochure are available in a range of fabrics (see below for Express Delivery range fabric).


Amongst our extensive range of fabrics we also offer a number of anti-microbial options that are effective against MRSA as well as E.coli, Listeria, Salmonella, Legionella and other similar bacteria.

These include


- Anti-microbial vinyl fabric
- Anti-microbial impervious fabric
- Intervene fabric, which is waterproof and has added stain repellent and anti-microbial properties



Anti-microbial vinyl fabric



Anti-microbial impervious fabric



Intervene fabric

For sample swatches or copies of the cleaning / care instructions on your chosen fabrics please contact our Sales Office on **01606 593456**.

Hygienilac for hygiene protection



All wooden frames are now supplied with the hygienilac coating. This coating incorporates an active agent which prevents bacteria surviving in contact with the coated surface.

Independent laboratory tests prove that hygienilac coated products have a 99% kill rate of these bacteria within a 24 hour period and hygienilac goes on killing bacteria for the life of the coated surface.


Express Delivery

EXPRESS DELIVERY


Certain items are available for immediate delivery from stock, and can be delivered within 15 days from receipt of order.

For these Express Delivery items look out for the special icon which is shown next to the relevant products.


Express seating products are supplied in the following vinyls:




Voyage 198 - Hyacinth



Voyage 105 - Indigo



Voyage 155 - Duck Egg



Berkshire Burghfield



Berkshire Newbury

Patient bedside furniture

Our high quality and practical patient furniture, designed for use in a ward setting, has been carefully constructed with unique options and to meet the needs of today's healthcare environments.



Patient bedside furniture

Enterprise tables12 - 13
Enterprise lockers14 - 15
600 range lockers16
A range lockers17
100 range combination lockers18
Elite patient chair range19 - 21
Baroque chair range22
Balero chair range23
Visitor seating24

Enterprise tables

The Enterprise™ tables are lightweight, durable and easy to clean. For complete flexibility they have been designed for the patient to use both in bed and when seated in a chair.



Colour options



Beech



Blue



Grey

EXPRESS DELIVERY

PT10K Enterprise overbed / chair table (beech)
MBC Plastic base cover



EXPRESS DELIVERY

Colour options



Beech





Blue



Grey

PT11K Enterprise split / tilting overbed / chair table (beech, blue or grey)*

Furniture item		
Product code	PT10K	PT11K
Height range	710mm - 1080mm	710mm - 1080mm
Table top	900mm x 420mm	900mm x 420mm



50mm twin castor wheels

Features:

- 50mm twin wheel castors for easy manoeuvrability
- Gas spring assisted height adjustment
- Excellent height range of 710 - 1080mm
- PT11 - Tilting top adjusted to a variety of angles using lightweight, easy to use ratchet system
- Smooth, non-absorbent surfaces for ease of cleaning
- Table top 900mm x 420mm
- Safe working load (evenly distributed) 25kg
- Product weight 17kg



Enterprise lockers

The Enterprise™ lockers are the first in the market to comprise of a rigid cast aluminium base and corner posts. This makes them extremely strong and durable, as well as giving them a novel contemporary appearance.

Features:

- One piece plastic moulded drawers and smooth surfaces for easy cleaning
- Bottle holder (holds 3 x 1.5L bottles)
- Free 65mm castors (2) and kick stop castors (2)
- Dual sided bag hooks and retractable towel rail
- Lightweight, easy to manoeuvre with ergonomic easily accessible drawers
- Sturdy corner buffers
- Large bottom drawer with internal shelf



Furniture item		
Product code	PL2	PL3
Height	850mm	850mm
Width	450mm	450mm
Depth	470mm	470mm



PL3 Enterprise 3 drawer patient bedside locker (beech)



PL3 Enterprise 3 drawer patient bedside locker (grey)

EXPRESS DELIVERY



PL2 Enterprise 2 drawer patient bedside locker (grey)



Plastic moulded drawers



PL2 Enterprise 2 drawer patient bedside locker (beech)



Removable drawers and internal shelf for easy cleaning

EXPRESS DELIVERY

Options:

- Lock (L-CAB)
- Master key (MK)
- Suited locking system (SL)
- Master key - suited locking system (MKSLS)
- Grand master key - suited locking system (GMSLS)
- Removable plastic top cover (PLRPC)

Add the relevant letter to the product code when ordering

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

600 range lockers

Flexible design, made to a robust construction with a smooth laminate finish including a moulded top.

Features:

- Moulded top
- 80mm 360 degree swivel castors and buffers
- Post formed smooth laminate sides
- Smooth laminate finish
- Towel rail to the rear
- Bowl holder optional
- Lock with two individual keys

Optional:

- Master keys
- Suited locking system
- Braked castors



DB602 Patient bedside locker / drop-down flap



DB652 Patient bedside robe locker / drop-down flap

Layout

- Self medicating compartment internal dimensions W: 416mm x D: 146mm x H: 220mm, drop down flap
- Through doors, accessible from the left and right
- Middle drawer
- 2 door cupboard with adjustable shelf
- Robe 255 wide with hanging rail

Furniture item		
Product code	DB602	DB652
Height	1050mm	1400mm
Width	455mm	710mm
Depth	455mm	515mm

A range lockers

Self medication lockers practically designed.

Layout

- Self medicating lockable compartment, drop-down flap, heavy duty hinge (A2, A3)
- Deep lockable top drawer (A4)
- Through doors, accessible from the left and right (A2 & A3)
- Middle drawer (A4)
- 2 door cupboard with adjustable shelf (All)
- Robe 255 wide with hanging rail (A3)
- Self medicating compartment internal dimensions W: 416mm x D: 146mm x H: 220mm
- Self medicating drawer dimensions W: 360mm x D: 330mm x H:290mm



TYPE A3 Self medication combination locker



TYPE A4 Self medication bedside locker / drawer

Furniture item			
Product code	TYPE A2	TYPE A3	TYPE A4
Height	975mm	1400mm	975mm
Width	460mm	715mm	460mm
Depth	515mm	515mm	515mm

EXPRESS DELIVERY

TYPE A2 Self medication bedside locker / drop-down flap

Features:

- 75mm 360 degree swivel castors
- Smooth laminate finish
- Towel rail to the rear
- Bowl holder
- Lock with two individual keys

Optional:

- Master keys
- Suited locking system
- Braked castors



To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

100 range combination lockers

High quality and versatile, covering various needs with many built in features.

Features:

- 75mm 360 degree swivel castors
- Smooth laminate finish
- Chrome handle
- Towel rail to the rear (bowl holder optional)

Optional:

- Braked castors

Layout

- 1 Drawer / 1 Door cupboard with adjustable shelf (DB103)
- 4 Drawers (DB101)
- Robe 255 wide with hanging rail (DB101, DB103)



Braked castor





Standard castor



DB101 Combination locker and robe



DB103 Combination locker

Furniture item		
Product code	DB103	DB101
Height	1400mm	1400mm
Width	715mm	715mm
Depth	515mm	515mm






Elite patient chair range

The Elite™ range has been designed using anthropometric data which highlights 13 critical dimensions to achieve a good sitting posture. These 13 dimensions are included in each of the 6 chair sizes that make up the Elite™ range. Sitting in the correct sized chair ensures the patient's comfort and wellbeing.



EXPRESS DELIVERY



Furniture item		
Product code	E600 with integral wing	E601 without wings
Seat height	482mm	482mm
Seat width	502mm	502mm
Seat depth	457mm	457mm
Seat height options:	16",17",18",19",20",21",22"	

PRESSURE REDUCING CUSHION AND HOUSEKEEPING WHEELS SUPPLIED AS STANDARD

E601 Elite chair without wings

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Elite patient chair range




E351 Elite medium chair with drop-down arms

Moving and handling

- The Elite chairs have an option of an adjustable height arm, aiding lateral transfer
- The housekeeping wheels allow the chair to be easily moved (when the chair is empty)
- Sacral gap allows easy access to fit hoist slings
- The arms allow easy grip with rounded ends and a stable base to push up from
- Splayed legs ensure greater stability and easy hoist access
- SWL 127kg (20 stone)

Tissue viability

- The visco-elastic seat cushion helps in the treatment and management of pressure ulcers with frequent monitoring and repositioning.
- Impressive pressure reduction testing
 - The 10° recline assists spinal alignment
 - Option of a seat board (SB) to allow specialist cushions to be used



Elite chair drop-down arm option

- Adjustable height arms aid lateral transfer
- Arms can be reduced to any point using an easy to use pull button

Optional

- Armpads
- Seat board



Seat board

Standard



Pressure reducing cushion



Housekeeping wheels



Size identification discs



00 Integral wing



01 Without wings



E200 Small Elite chair



E300 Medium Elite chair



E400 Large Elite chair

Furniture item			
Product code	E200	E300	E400
Seat height	432mm	451mm	470mm
Seat width	470mm	502mm	527mm
Seat depth	438mm	457mm	476mm

Infection control

The Elite chairs have been specifically designed to be easy to clean. Features include:

- Saddle stitched seams to reduce the risk of fluid ingress
- Removable seat cushion for optimal cleaning access
- Sacral gap to allow easy cleaning and access
- Optional anti-microbial fabrics
- Anti-microbial hygienilac coating on the wooden frame

SEE PAGE 26 FOR THE ELITE BARIATRIC X-LARGE CHAIR

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Baroque chair range

Various options provide a chair that satisfies every need.

Options:

- Pressure reducing seat cushion (P)
- Wings (W)
- Armpads (A)
- Housekeeping wheels (H)
- Filled-in sides (F)
- Seat height standard 19"
- Seat height optional 16", 17", 18" or 20"
- Medium and low back, see page 33

Add the relevant letter to the product code when ordering

Features:

- Sacral gap and removable seat cushion allow easy cleaning and access
- An optional pressure reducing seat cushion for the management and prevention of pressure ulcers
- Optional housekeeping wheels make the chair easy to move reducing moving and handling risks
- The seat cushion has two seat depth positions for increased patient comfort and well being
- All wooden frames are coated in hygienilac lacquer and anti-microbial fabrics are available



EXPRESS DELIVERY

BR300 Baroque high back chair



W Wings



F Filled-in sides




A Armpads
P Pressure reducing cushion



H Housekeeping wheels

Furniture item	
Product code	BR300
Height	1219mm
Width	610mm
Depth	686mm
Seat height	406 - 508mm
Seat width	483mm
Seat depth	432 or 483mm



To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Balero chair range

Attractive and practical design with a robust frame.



W Wings



A Armpads
P Pressure reducing cushion




H Housekeeping wheels



J Height adjusters

Furniture item	
Product code	BL300
Height	1219mm
Width	610mm
Depth	686mm
Seat height	406-508mm
Seat width	483mm
Seat depth	483mm



Features:

- Sacral gap and an optional removable seat cushion allow easy cleaning and access
- An optional pressure reducing seat cushion for the management and prevention of pressure ulcers
- Optional housekeeping wheels make the chair easy to move reducing moving and handling risks
- All wooden frames are coated in hygienilac lacquer and anti-microbial fabrics are available

EXPRESS DELIVERY

BL300 Balero high back chair



Options:




- Pressure reducing seat cushion (P)
- Wings (W)
- Armpads (A)
- Housekeeping wheels (H)
- Removable seat cushion
- Seat height standard 19"
- Seat height optional 16", 17", 18" or 20"
- Medium and low back, see page 34

Add the relevant letter to the product code when ordering

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Visitor seating

Varied options are available for visitors in the ward environment.

Furniture item			
Product code	F833	F835	VENUS
Height	850mm	850mm	835mm
Width	550mm	550mm	470mm
Depth	530mm	530mm	510mm
Seat height	483mm	483mm	440mm
Seat width	500mm	500mm	470mm
Seat depth	570mm	570mm	400mm



EXPRESS
DELIVERY

F833 Finland stacking chair



EXPRESS
DELIVERY

F835 Finland stacking chair with straight arms



VENUS Visitor chair

Bariatric / recliners & specialist seating

All our chairs are made to a high quality robust construction, with flexible additions available.

Bariatric / recliners & specialist seating

Bariatric chair range	26
Recliners	27
Specialist seating	28



Bariatric chair range




All our chairs are made to a high quality robust construction, with flexible additions available.

Bariatric features:

- Pressure reducing seat cushion helps in the treatment and management of pressure ulcers with frequent monitoring and repositioning
- Housekeeping wheels for easy manouvability when empty
- Sacral gap allows easy access to fit hoist slings
- Sturdy and robust frame
- Removable seat cushion
- SWL 254kg (40 stone)

Elite features:

- Pressure reducing seat cushion as standard
- Housekeeping wheels as standard
- Sacral gap allows easy access to fit hoist slings
- The arms allow easy grip with rounded ends and stable base to push up from
- The Elite chairs have an option of an adjustable arm, aiding lateral transfer
- SWL 159kg (25 stone)

Furniture item			
Product code	E525	E535	AC300HW
Seat height	495mm	495mm	508mm
Seat width	572mm	572mm	711mm
Seat depth	508mm	508mm	610mm



AC300HW Bariatric chair & housekeeping wheels



E525 X-Large Elite chair






E535 X-Large Elite chair with drop-down arms

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Recliners

An extremely comfortable recliner available in a manual or electric version.

							
Furniture item	Manual		Electric		Bariatric		
Product code	RR2	RR3	RR1	RS1	RS2	DR1	RB1
SWL (kg)	121	159	121	121	159	254	318
SWL (stone)	19	25	19	19	25	40	50
Motor	n/a	n/a	1	2	3	3	3
Seat height	457mm	457mm	457mm	457mm	457mm	457mm	457mm
Seat width	508mm	559mm	508mm	508mm	559mm	711mm	711mm
Seat depth	508mm	508mm	508mm	508mm	508mm	508mm	508mm



RR1 Renray electric recliner



Cushion back

RR2 Renray manual recliner



Waterfall back

The waterfall back available on the Renray recliner is fibre filled and can be adjusted to suit the user's needs. This allows for excellent support in any position, whether seated or reclined.

Features:


- Braked castors (BC)
- Pressure reducing cushion (P)
- Back height 762
- The RS1 is a 2 motor 4th position chair with the 4th position being Trendelenburg
- The RR2, RR3, RR1, RS1 and RS2 have a waterfall style back
- The DR1 and RB1 have a cushion style back

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Specialist seating

The Delta chair with gas spring adjustment, enables fluid precise movement of the tilt-in-space facility by the release of a lever situated on the push handle. The tilt-in-space gently moves back and forth as directed by the carer with minimal effort whilst the leg rest can be elevated by a separate release mechanism at the front on the chair.



Furniture item	
Product code	DL002A
Seat depth	380-480mm
Seat width	450mm
Seat height	530mm

DL002A Delta 450 with gas adjustment

Features:

- Pressure relieving Reflexion™ foam
- Multi-positioning
- Seat length adjustment
- Arcing leg rest motion without compromising footboard positioning
- Safe working load 180kg / 28 stones
- Optional 550mm seat width DL004A

Waiting area

These items are designed to complement the rest of our range to ensure that the overall look is consistent and there are appropriate pieces to fill every need. All have been made by our craftsmen with both style and durability in mind.






Waiting area

Finland chairs (stackable)30
Venetian chairs (stackable)31
Venetian chairs32
Baroque chairs33
Balero chairs34
Coffee tables35

Finland chairs (stackable)

A modern style available in a choice of contemporary colours.

Furniture item			
Product code	F833	F834	F835
Height	850mm	850mm	850mm
Width	550mm	550mm	550mm
Depth	530mm	530mm	530mm
Seat height	483mm	483mm	483mm
Seat width	500mm	500mm	500mm
Seat depth	570mm	570mm	570mm

OPTIONAL LINKING MECHANISM



EXPRESS DELIVERY

F833 Finland stacking chair



EXPRESS DELIVERY

F834 Finland stacking chair with arms



EXPRESS DELIVERY

F835 Finland stacking chair with arms





N.B. Recommended stacking six chairs maximum



To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Venetian chairs (stackable)

Comfortable and stylish with a robust frame, suitable for most areas.

Furniture item		
Product code	V010	V020
Height	820mm	820mm
Width	580mm	580mm
Depth	550mm	550mm
Seat height	457mm	457mm
Seat width	495mm	495mm
Seat depth	482mm	482mm

V020 HAS A CURVED BACK TO ENSURE MAXIMUM COMFORT WHEN SEATED FOR LONG PERIODS



V010 Venetian side chair



V020 Venetian side chair with arms



N.B. Recommended stacking three chairs maximum



To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Venetian chairs





Comfortable and stylish with a robust frame, suitable for most areas.

Features:

- Sacral gap on all chairs for easy cleaning
- Strong and durable laminate frame
- V230 Venetian chairs have easy to grip arms to aid rising
- V300 and V080 stackable 3 high



V300 Venetian low back chair with arms

Furniture item				
Product code	V230	V300	V080	V500
Height	1194mm	864mm	864mm	-
Width	610mm	584mm	508mm	-
Depth	660mm	610mm	610mm	-
Seat height	457mm	457mm	457mm	457mm
Seat width	553mm	483mm	483mm	-
Seat depth	457mm	432mm	483mm	-



V230 Venetian high back chair



V080 Venetian low back chair without arms



V500 Venetian footstool

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Baroque chairs

Practical and stylish design in a range of options.






BR105 Low back chair with arms



BR140 Medium back chair with arms



BR100 Low back chair without arms

Furniture item			
Product code	BR140	BR100	BR105
Height	1041mm	940mm	940mm
Width	610mm	533mm	610mm
Depth	660mm	635mm	635mm
Seat height	406-508mm	406-508mm	406-508mm
Seat width	483mm	483mm	483mm
Seat depth	432 & 483mm	432 & 483mm	432 & 483mm

Features:

- Sacral gap and removable seat cushion allow easy cleaning and access
- The seat cushion has two seat depth positions for increased patient comfort and well being
- All wooden frames are coated in hygienilac lacquer and anti-microbial fabrics are available
- Seat height standard 19"
- Seat height optional 16", 17", 18" or 20"

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Balero chairs



Comfort and support with a clean and practical design.

Features:

- Sacral gap allows easy cleaning
- All wooden frames are coated in hygienilac lacquer and anti-microbial fabrics are available
- Seat height standard 19"
- Seat height optional 16", 17", 18" or 20"



BL105 Balero low back chair with arms

Furniture item		
Product code	BL105	BL140
Height	940mm	1041mm
Width	610mm	610mm
Depth	635mm	635mm
Seat height	406-508mm	406-508mm
Seat width	483mm	483mm
Seat depth	483mm	483mm



BL105 Balero low back chair with arms








BL140 Balero medium back chair

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Coffee tables

Versatile and stylish, available in an attractive smooth laminate finish.

Furniture item					
Product code	CTC	CTCL	CTR	CTVS	CTU
Height	640mm	457mm	460mm	440mm	520mm
Width	600mm DIA	600mm DIA	920mm	560mm	600mm
Depth	-	-	460mm	1000mm	600mm



CTC Circular coffee table



CTCL Circular coffee table low height



CTR Rectangular coffee table



CTVS Oval coffee table with shelf



CTU Uni coffee table



- Colour options:
- Standard:
- Beech
- Optional
- Maple
 - Pippy Oak
 - Walnut
 - Cherry

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Dining furniture

The items in the comprehensive dining range come in various shapes to ensure whatever the need, it can be met. Our dining chairs offer flexible options depending upon the different requirements of the end user.



Dining furniture

Bistro range38-39
Chelford dining chairs40
Dining tables41

Bistro range

Modern and stylish looks combine with practical comfort.

Features:

- Beech finish
- Seating available with or without seat pads
- Chrome frames



MARSWAWP Mars bistro chair with seat pad










MARS Mars bistro chair



MARSWA Mars bistro chair with arms



MARSHSWP Mars bistro high stool with seat pad

Furniture item							
Product code	MARS	MARSWAWP	MARSHSWP	LUNARVARPED	LUNARVAR	LUNARVAS	TRUMPETB2
Height	840mm	810mm	1136mm	710mm	710mm	710mm	720mm
Width	500mm	470mm	544mm	760mm	700mm DIA	700mm	700mm DIA
Depth	530mm	470mm	635mm	1220mm	-	700mm	-
Seat height	460mm	460mm	759mm	-	-	-	-
Seat width	450mm	450mm	450mm	-	-	-	-
Seat depth	457mm	457mm	457mm	-	-	-	-

AVAILABLE WITH A
SMALL SQUARE OR CIRCULAR TOP

LUNARVARPED Rectangular bistro table



TRUMPETB2 Circular bistro table




To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

40

Chelford dining chairs

Contemporary and versatile with a variety of options.



Furniture item			
Product code	CHE	CHEWA	CHEWAS
Height	889mm	889mm	889mm
Width	495mm	495mm	495mm
Depth	470mm	470mm	470mm
Seat height	457mm	457mm	457mm
Seat width	457mm	457mm	457mm
Seat depth	432mm	432mm	432mm



CHE Chelford chair



CHEWA Chelford chair with arms






CHEWAS Chelford chair with arms & skids

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Dining tables

A variety of stylish shapes and designs to cover all dining requirements, in a smooth, heat resistant laminate finish.

Furniture item			
Product code	DTC	DTS	DTR
Height	760mm	760mm	760mm
Width	1070mm DIA	910mm	1520mm
Depth	–	910mm	910mm



DTC Circular dining table



DTS Square dining table



DTR Rectangular dining table

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

- Colour options:
- Standard:
- Beech
- Optional
- Maple
 - Pippy Oak
 - Walnut
 - Cherry

Sofas, upholstered chairs & tub style seating

Our range of seating has been carefully selected to withstand day to day demands. Every piece of furniture in this range is constructed using strong, durable materials with many featuring built in extras.






Sofas, upholstered chairs & tub style seating

Belton range44
Harvard range45
Milan range46
Kingsley range47
Chelford range48
Oulton range49

Belton range

Stylish looks and modern practicality, in a high quality leather effect finish.

Furniture item			
Product code	BEL	BEL2S	BEL3S
Height	940mm	940mm	940mm
Width	910mm	1470mm	2010mm
Depth	970mm	970mm	970mm
Seat height	480mm	480mm	480mm
Seat width	560mm	1120mm	1650mm
Seat depth	530mm	530mm	530mm



BEL Belton chair






BEL2S Belton 2 seater sofa

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Harvard range

Timeless elegance with a practical and homely feel.

Furniture item			
Product code	HAC	HA2S	HA3S
Height	915mm	915mm	915mm
Width	940mm	1423mm	1918mm
Depth	889mm	889mm	889mm
Seat height	432mm	432mm	432mm
Seat width	483mm	940mm	1423mm
Seat depth	559mm	559mm	559mm



HAC Harvard chair



HA2S Harvard 2 seater sofa



To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Milan range



Attractive seating with flexible and practical adaptability.



MIL2SB Milan 2 seater sofa bed






MILCHB Milan chair bed

Furniture item		
Product code	MILCHB	MIL2SB
Height	800mm	800mm
Width	965mm	1562mm
Depth/bed length	737mm	1829mm
Seat height	419mm	419mm
Seat width	750mm	1500mm
Seat depth	495mm	495mm

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Kingsley range

Sophisticated and simple shapes help create a stylish look.

Furniture item			
Product code	KINH	KIN	KIN2S
Height	1041mm	838mm	838mm
Width	610mm	610mm	1219mm
Depth	711mm	711mm	711mm
Seat height	460mm	460mm	460mm
Seat width	508mm	508mm	1041mm
Seat depth	482mm	482mm	482mm



KINH Kingsley high back chair



KIN2S Kingsley low back 2 seater sofa






KIN Kingsley low back chair

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Chelford range

Flexible options come in a contemporary shape.

Furniture item			
Product code	CHEL	CHELL	CHELT
Height	1060mm	890mm	890mm
Width	780mm	780mm	780mm
Depth	720mm	720mm	720mm
Seat height	508mm	508mm	432mm
Seat width	520mm	520mm	520mm
Seat depth	500mm	500mm	500mm



CHEL Chelford chair



CHELL Chelford low back chair





CHELT Chelford tub chair

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Oulton range

Modern designs help to create a stylish look.

Furniture item		
Product code	OUL	OULL
Height	1100mm	810mm
Width	720mm	740mm
Depth	760mm	670mm
Seat height	510mm	510mm
Seat width	540mm	540mm
Seat depth	510mm	510mm



OUL Oulton high back chair



OULL Oulton low back chair

To make an enquiry or to purchase this product call 01606 593456 or email sales@renrayhealthcare.com

Your next step

Contact us

For further information, including product availability, the latest colour options and prices, and to order any of the products in this brochure; or if you have any queries about other Renray products and services, contact our Sales Office on:

Tel: 01606 593456

Email: sales@renrayhealthcare.com

Fax: 01606 861354

Our experienced team will be pleased to discuss your requirements. Also our nationally based Area Sales Managers are always available to visit you for a no obligation discussion and to help advise you, should you prefer.

Renray, furniture that's the perfect fit for any environment

If you would like further information about other products and services available from the Renray Group, please contact us via the details above.

Challenging Environments

Care Homes

Dementia

Express Delivery

Hotel & Leisure

Housing Associations

Your next step

Contact us51